


OCA CORES: the Optimal Solution for Small, Osteochondral Defects

JRF Ortho offers **Precut Fresh Osteochondral Allograft (OCA) Cores** for the treatment of full thickness osteochondral lesions.

Available in 10mm, 12mm and 16mm diameter sizes, The Cores provide complete replacement with viable intact osteochondral allograft without donor site morbidity.

PRECUT FRESH OCA CORES:
A GREAT ALTERNATIVE FOR
FOCAL LESIONS OF 20 MM OR LESS



Indications

- ▶ Knee
- ▶ Shoulder
- ▶ Foot / Ankle
- ▶ Elbow

OPTIMAL AND CONVENIENT TREATMENT

Structurally intact: contains a mature, viable, hyaline cartilage matrix with bone for superior biomechanical and durable repair

Single step treatment with no patient size matching required

Proprietary processing methods shown to maintain 80% absolute cartilage viability¹

Utilizes existing instrumentation – disposable Arthrex® OATS® sets*

* OATS® is a registered trademark of Arthrex, Inc
Disposal set ordered separately through Arthrex, Inc
Data on file, JRF Ortho

10 MM AND 12 MM CORE IMPLANTATION

technique guide

PATIENT PREPARATION

Follow the standard Autograft OATS® technique (Step 7–9 in the Autograft OATS technique guide) for the preparation of the recipient implantation site for the Precut Fresh OCA Cores.

GRAFT PREPARATION

The following are alternative steps required for successful implantation of the Precut Fresh OCA Cores:

- ▶ Prepare the length of the allograft core to match the depth of the harvested site
- ▶ Load The Core in clear plastic sleeve
- ▶ Face the cartilage side of the Fresh OCA Core up or away from beveled end of the plastic sleeve
- ▶ With the core extruder knob handle fully retracted, slide the plastic sleeve with The Core onto the blue handled donor harvester device widest side first
- ▶ To insert The Core into the donor harvester, ensure The Core lines up with the chamber of the blue handled harvester device and then gently seat The Core into the harvester
- ▶ Once The Core is securely inserted into harvester device, place the beveled end of the plastic sleeve in the recipient site and deploy white extruder knob to fully seat The Core into place

JRF ORTHO PRECUT FRESH OCA CORES

| ALLOGRAFT | VOLUME | PART CODE |
|---------------------------|--------|-----------|
| Precut Osteochondral Core | 10 mm | 45647010 |
| Precut Osteochondral Core | 12 mm | 45647012 |
| Precut Osteochondral Core | 16 mm | 45647016 |

ARTHREX® INSTRUMENTATION

| INSTRUMENT | SIZE | ID |
|---------------------|-------|--------------|
| Single-Use OATS Set | 10 mm | ABS-8981-10S |
| Single-Use OATS Set | 12 mm | ABS-8981-12S |
| Single-Use OATS Set | 16 mm | ABS-1981-16S |

**16
MM**

16 mm core implantation: Follow the technique for the Arthrex Single-Use OATS Set, 16 mm ABS-1981-16S

For instrument info: Arthrex.com

© JRF ORTHO 2023

TO ORDER

JRFORTHO.org/order
877-255-6727



7245 S Havana St, Ste #300
Centennial, CO 80112
M 877-255-6727 F 303-649-1658
JRFORTHO.org

